HYSIS MEDICAL

LAVE-BASSIN LPC

BED-PAN WASHER LPC

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NEW MODEL:

- . THE BED-PAN CAN ONLY BE EMPTIED WITH THE LID CLOSED (ANTI-SPLASH)
- . BED-PAN TILTING DURING WASHING
- . COMPACT SIZE
- . ADAPTABLE TO ANY ITEM TO BE WASHED

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LPC bed-pan washer wall mounted with sliding door shut

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LPC bed-pan washer floor mounted with sliding door open

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● AUTOMATIC BED-PAN WASHER

Thermal disinfection
Automatic opening LPC
Foot operated in option

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● The automatic bed-pan washers LPC allow to empty, wash with high pressure, rinse and disinfect all type of basins, urinals and bottles. Totally built in stainless steel AISI 304. The sink is made of AISI 316.

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Easy to use

- Front loading with:
 - fix support allowing all models of basins or square commode containers
 - 1 or 2 independents nozzles for urinals or bottles
 - 1 cover handling device
 - accessories upon request
- Cycle starting operated through the front panel or by foot action (LCA model) The front control panel gives clear visual indications about:
 - cycle operation
 - temperature of disinfection
 - occasional dysfunction
 - number of cycles

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Easy to handle and evolutive

Microprocessor cycle management with user defined possibilities upon request:

- by prompt activation of the menu on front panel (Double washing, double rinsing, etc)
- software modification (cycle duration, automatic opening at the end of the cycle, etc)

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High technical liability

- Technical components hardly tested
- Easy maintenance thanks to the front and lateral panels (Access to the main hardware components)
- Water softening device
- Electro valves with stainless steel protecting grill

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Safe

- Cycle registration on paper thanks to a PC connection
- Lockable cover
- Emergency stop and reset

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Noiseless

Sound proof allowing a noise at maximum level 50 dB (A).

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● WASHING UNIT:

The bed-pan washers can be associated with emptiers, sinks, storage platforms with sliding doors and/or brackets. The hydraulic/electric connections between the various components are preassembled in the factory and do not require any particular work during installation.

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● EMPTIER:

Totally made of AISI 304 stainless steel with 400x400 mm pressed sink fitted with rinsing nozzles controlled by a quick-flow tap. Fold-down support grid for buckets and pots. Complete with tap for filling them. Opening panel on the front with small bay available and drain coupling.

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● WC VERSION:

Particularly suitable for restructuring existing rooms.

Version with:

- stainless steel WC
- ceramic WC
- complete with storage tank and WC flushing tank or, alternatively, with a tank with external controls.

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FOR FURTHER INFORMATION SEE THE BROCHURE ABOUT THE 302 PROGRAMMABLE ELECTRONIC CONTROL UNIT

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Serial port for unloading the data stored on personal computer and printer (access for maintenance)

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


● Evolved control of the temperature and disinfection times.

Hatch opening/closing system with electric lock that prevents it being opened during the cycle and microswitch that prevents enabling of the cycle when the hatch is open. Friction damper that keeps the hatch stable in the open position.

Audible alarm in the case of faults during the cycle.

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● *Washing/disinfection system: the slipper bedpan is washed by optimising the action of the nozzles with a conical jet and cutting washing action as the bedpan is tilted. This gives considerable energy saving (water consumption over the same time period is approximately 50% lower than with rotary action jets). The disinfection cycle is done by output of flowing steam. The temperature read in the chamber can be regulated between 80 and 95°C and maintained for an adjustable time. The addition of bleaching agent to the rinsing water is available as an option. Dilution is adjustable.*

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● *The thermal disinfection washer with self-generated steam is normally fitted with a liquid water-softening device with dosing pump and level gauge.*

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


● AVAILABLE MODELS WITH REDUCED POWER


Lave-bassin LPC mural
LPC bed-pan washer wall mounted

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● Electrical and hydraulic connections

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

● Certification

- The washer is manufactured and tested according to the CEI-IEC standards in force and bears CE marking
- Disinfection has been tested at various microbiological laboratories (certification available on request)

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Hysis Medical consistently with its policy of continuous improvement on its own products, reserves the right to make changes at any time without prior notice. Therefore, all related data must be regarded as approximate and unbinding.